

Cambiador de muestras rotatorio TW 7200

Cambiadores de muestras



Imagen referencial de la serie (fabricante).

Variantes	10 (ver tabla de variantes)
Para Chile	Fuente 115-240 V, 50/60 Hz, 30 VA: válida para la red chilena (220 V / 50 Hz)
Interfaces	2 x RS-232-C, 1 x I/O, 2 x pumps, 1 x stirrer, 1 USB-A (Update only), 1 USB-Mini-B (PC), 1 x Ethernet
Dimensiones	Basic unit: 143 x 620 x 505 mm (W x H x D)Sample tray: 450 mm
Peso	10.3 kg Examples with sample trays:Sample tray with 16 beaker glasses 150 ml: approx.2.6 kgSample tray with 24 beaker glasses 50 ml: approx 2.4 kgSample tray with 24 rea...
Alimentación	115 - 240 V~, 50/60 Hz, 30 VA
Agitador magnético	2 fixed magnetic stirrers

Especificaciones

Interfaces	2 x RS-232-C, 1 x I/O, 2 x pumps, 1 x stirrer, 1 USB-A (Update only), 1 USB-Mini-B (PC), 1 x Ethernet
Dimensiones	Basic unit: 143 x 620 x 505 mm (W x H x D)Sample tray: 450 mm
Peso	10.3 kg Examples with sample trays:Sample tray with 16 beaker glasses 150 ml: approx.2.6 kgSample tray with 24 beaker glasses 50 ml: approx 2.4 kgSample tray with 24 rea...
Alimentación	115 - 240 V~, 50/60 Hz, 30 VA
Agitador magnético	2 fixed magnetic stirrers
Motor horizontal	Horizontal motor for driving a second sample beaker row
Gradillas	From 12 - 48 positions
Normas	EN 61010-1 (seguridad)

Variantes

Código	Variante	Muestras (gradilla)	Bomba de enjuague	Contenido
285228000	TW 7200	—	—	Basic unit TW 7200 with two integrated magnetic stirrer, connection cable for rod stirrer, USB connection cable TZ 3887, power supply 100- 240 V
285228010	TW 7200-16	16	—	Basic unit with sample rack TZ 4016 for 16 samples, without beakers, incl. titration head TZ 1463 and power supply 100- 240 V. Suitable beakers are 150 ml TZ 1784, TZ 1788 or 250 ml TZ 1787, TZ 1780
285228020	TW 7200-18	18	—	Basic unit with sample rack TZ 4018 for 18 samples, without beakers, incl. titration head TZ 1464 and power supply 100- 240 V. Suitable beakers are TZ 1765 (100 ml tall form without spout)
285228030	TW 7200-24	24	—	Basic unit with sample rack TZ 4024 for 24 samples, without beakers, incl. titration head TZ 1469, power supply 100- 240 V. Suitable sample beakers/vessels are 50 ml TZ 1783 / 75 ml TZ 1786

Código	Variante	Muestras (gradilla)	Bomba de enjuague	Contenido
285228040	TW 7200-COD	—	—	Basic unit with sample rack TZ 4040 for 24 CSB-samples according to DIN 38 409, without vessels, incl. Titration head TZ 1461, redox electrode Pt 5901, rod stirrer TZ 1846, titration tip TZ 1648, power supply 100- 240 V
285228050	TW 7200-16-MP	16	MP 25	Basic unit with sample rack TZ 4016 for 16 samples, without beakers, incl. titration head TZ 1467, rinsing pump MP 25, power supply 100- 240 V. Suitable beakers are 150 ml TZ 1784, TZ 1788 or 250 ml TZ 1787, TZ 1780
285228200	TW 7200-25	25	—	Basic unit with sample rack TZ 4025 for 25 samples, without beakers, incl. titration head TZ 1463 and power supply 100- 240 V. Suitable beakers are 150 ml TZ 1784, TZ 1788 or 250 ml TZ 1787, TZ 1780

Código	Variante	Muestras (gradilla)	Bomba de enjuague	Contenido
285228210	TW 7200-25-MP	25	MP 25	Basic unit with sample rack TZ 4025 for 25 samples, without beakers, incl. titration head TZ 1467, rinsing pump MP 25, power supply 100-240 V. Suitable beakers are 150 ml TZ 1784, TZ 1788 or 250 ml TZ 1787, TZ 1780
285228220	TW 7200-28	28	—	Basic unit with sample rack TZ 4028 for 28 samples, without beakers, incl. titration head TZ 1464 and power supply 100- 240 V. Suitable beakers are TZ 1765 (100 ml tall form without spout)
285228230	TW 7200-38	38	—	Basic unit with sample rack TZ 4038 for 38 samples, without beakers, incl. titration head TZ 1469, power supply 100- 240 V. Suitable sample beakers/vessels are 50 ml TZ 1783 / 75 ml TZ 1786

Ficha preparada por Cromtek con información del fabricante (xylemanalytics.com (web oficial SI Analytics), 2026-10-08). Datos sujetos a confirmación; un ejecutivo de Cromtek le enviará precio y plazo de entrega.

Cromtek · www.cromtek.cl · contacto@cromtek.cl · (+56 2) 3378 0171 · +56 9 4413 5160